



Development of a multi-parameter flow cytometric assay for the detection of primed T cells

A.D. Kelleher

**Head Immunovirology and Pathogenesis
program**

**National Centre HIV Epidemiology and
Clinical Research**





Elispot/ ICC

Short term assays: detect immediate cytokine production
good at detecting effector cells but perhaps not central memory cells

Proliferation assays

ability of cells to proliferate, but no insight into function or phenotype
CFSE labeling allows delineation of cell type proliferating



Objectives



- **Develop assay that allows detection of small numbers of primed or quiescent cells**
- **Assay should provide insight into function and phenotype**
- **It should be relatively simple to set up with relatively little *in vitro* manipulation**








Why?



SIV Gag-specific Response by IFN γ ELISPOT




Control of pathogenic SHIV challenge by DNA/FPV prime/boost immunisation





Week following SHIV_{mn229} IR challenge



Survival of monkeys post challenge








An assay to unmask priming

- **Combination of 2 Robust methods:**
 - CFSE staining to show proliferation
 - ICC assay to detect Ag specific responses (cytokine)
- **Addition of cytokine to boost small responses**
 - **IL-2, IL-7, IL-15:**
 - biological characteristics somewhat similar
 - act as T cell growth and survival factors
 - share the common cytokine receptor gamma chain (CD132)
 - implicated in memory cell generation
 - produced by different cells
- **c/w** Amplispot Nixon lab Effector cells
Cultured Elispot: Hill lab Non amplified



Basic proposed methodology



Parameters to be considered in initial work-up



•With and Without Cytokine:

- IL-2, IL-7, IL-15,

•Varying Concentrations of Cytokine

- (IL-2, IL-7, IL-15 & IL-7+ IL-15):

•Combinations of cytokine

•Length of incubation

•Varying time of cytokine addition:

- day 0 vs day 3 vs day 5

•Gating strategies

Other issues

- Optimal antigen

CMV-lysate/CEF

Measles



Vaccinia

- Frozen PBMCs

- 24 vs 96 well plate

- Retrieval of cells post EDTA











Assay Gating

CD4+

CD8+



IFN- γ

**CD4+CFSE
dim**

**CD4+CFSE
bright**

IL-2

**CD8+CFSE
dim**

**CD8+CFSE
bright**





Comparison of CD4 & CD8 responses to using 10 ng/ml IL-7 added at Day 0





No
Antigen

CD4+



CD8+



Vero



Measles



Non-specific background problematical







CD4 responses at varying concentrations of IL-7 added at Day 0

Vero



Measles

IL-7



10ng/ml



1ng/ml




0.1ng/ml





Non-specific proliferation is increased when cytokine is added at day 3 c/w day 5

↑
% Proliferation (CFSE dim)



Culture Conditions





DNA prime/rFPV boost candidate vaccine trial

Immunisation dose and schedule

DNA; 1mg (IMI) on day 0 and week 4

rFPV; 5×10^7 pfu/mL (IMI) on week 8

Immunogenicity By ELISpot /ICC/Proliferation at wks 0, 8, 9, 12, 24


CAMPICC assay: conducted week 8 after 2 DNA vaccines





DNA prime/rFPV boost candidate vaccine trial

Mean change from baseline ELIspot Response to Gag Pool 122






DNA prime/rFPV boost candidate vaccine trial





Lymphoproliferative responses

Mean change positive p24 SI





Vaccine responses appear to be CD4+ T cell mediated



Weeks of trial






CAMPICC assay: CD4+ T cell responses:

Vaccine recipient: #125, Wk8 Post DNA x2

Gag

Gag + IL-15





CAMPICC assay: CD4+ T cell responses:

Placebo recipient: #115, Wk8 Post DNA x2


Gag

Gag + IL-15



Day 10 + 0 + Gag

+ 0 + Gag



IL-2

#125 - Vaccine Recipient - Responder

Week 8 (post DNA x2): ICC Results

Gag+(CD28+CD49d)





Background values subtracted



CD4

CD8

CD-69



IFN-gamma



IL-2




#125 - Vaccine Recipient - Responder Week 9 (Post FPV Boost): ICC Results Gag+(CD28+CD49d)




Background values subtracted

CD-69



CD4



CD8



IFN-gamma



IL-2






CAMPICC assay: CD4+ T cell responses:

Vaccine recipient: #125, Wk8 Post DNA x2

Gag

Gag + IL-15



#115 - Vaccine Recipient - Placebo recipient Week 8(Post DNA x2): ICC Results







Gag+(CD28+CD49d)

Background values subtracted

CD4

CD8

CD-69



#115 - Vaccine Recipient - Placebo recipient

Week 9 (Post FPV Boost): ICC Results



Gag+(CD28+CD49d)





CD4

CD8

CD-69



IFN-gamma



IL-2





CAMPICC assay: CD4+ T cell responses:

Placebo recipient: #115, Wk8 Post DNA x2


Gag

Gag + IL-15



Day 10 + 0 + Gag

+ 0 + Gag




IL-2



VACCINE vs PLACEBO

GAG culture re-challenged with GAG

% CD4 proliferating and producing...



Placebo n=4

Vaccine n=11






VACCINE vs PLACEBO


GAG culture + IL-15 re-challenged with GAG

% CD4 proliferating and producing...


IFN- γ



IFN- γ & IL-2



IL-2



- + 109
- ✕ 110
- ▲ 111
- 113
- ▼ 114
- 115
- ◆ 116
- 117
- ◇ 118
- △ 119
- 120
- 121
- ⊗ 123
- ▽ 124
- 125
- ◇ 127

Placebo n=4

Vaccine n=11






VACCINE vs PLACEBO

GAG culture + IL-15 re-challenged with GAG

% CD8 proliferating and producing...





CD4 Responses

Week 8 & 9 ICC results vs CAMPICC assay



Vaccine

Placebo

IFN- γ



IL-2





CAMPICC assay: Results

Performed on last 16 individuals recruited

Evaluable data on 15: 11 active, 4 placebo

4 recipients of active vaccine had either IL-2 or IFN- γ production from CD4 that was $>$ mean +2 SD of the placebo responses,

1 recipient had CD8 response that fulfilled the same criteria

DNA may be priming 1/3 of the recipients

Not all of these boosted with rFPV

Suggests both components sub-optimal



Confirm by conducting optimised assay on all 24







ICS after two DNA vaccines: Macaques

IFN- γ



CD4+



CD8+







ICS following 10 day incubation with CFSE and IL7 (1ng/mL): macaque

IFN- γ

CD3+ CD4+ IFN γ +
Proliferating cells



CD3+ CD8+ IFN γ +
Proliferating cells





Comments/Conclusions

Assay initially worked up on robust responses:

? Appropriate methodology for borderline responses?

Unable to be validated in trial already conducted because of ineffective *in vivo* boost

Responses dominated by CD4, surprising, esp since using IL-15 (CD8 supporting cytokine) but is consistent with ICC data and data from similar vaccine approaches

IL-15 most effective cytokine so far in humans ? in macaques

Exploring other cytokines and cytokine combinations with Ag that induce less robust responses

Repeat analysis with fully optimized assay on cryo-preserved cells from this trial and proposed high dose phase I trial

Eliminate peptides driving “non-specific”/ cross reactive responses





Acknowledgements

IVPP/NCHECR

Mee Ling Munier

Susanna Ip

John Zaunders

David Cooper

University of Melbourne

Jane Dale

Rob de Rose

Stephen Kent

**This work was partially funded by NIH HVDDT award N01-AI-05395
and NHMRC Program Grant**

**NCHECR is funded by the Australian Commonwealth Department
of Health and Ageing**

